



Recombinant Human Tripartite motif-containing protein 65 (TRIM65)

Product Code	CSB-EP750884HU-B
Abbreviation	TRIM65
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PJ69
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAQLLEEKL TCAICLGLYQ DPVTLPCGHN FCGACIRDWW DRCGKACPEC REPFDPGAEL RRNVALSGVL EVVRAGPARD PGPDPGPGPD PAARCPRHGR PLELFCRTEG RCVCSVCTVR ECRLHERALL DAERLKREAO LRASLEVTQQ QATQAEGQLL ELRKQSSQIQ NSACILASWV SGKFSSLLQA LEIQHTTALR SIEVAKTQAL AQARDEEQRL RVHLEAVARH GCRIRELLEQ VDEQTFMQES QLLQPPGPLG PLTPLQWDED QQLGDLKQLL SRLCGLLLEE GSHPGAPAKP VDLAPVEAPG PLAPVPSTVC PLRRKLWQNY RNLTDFPVSA NRHFYLSRQD QQVKHCRQSR GPGGPGSFEL WQVQCAQSFQ AGHHYWEVRA SDHSVTLGVS YPQLPRCRLG PHTDNIGRGP CSWGGLCVQED SLQAWHNGEA QRLPGVSGRL LGMDLDLASG CLTFYSLEPQ TQPLYTFHAL FNQPLTPVFW LLEGRTLTLG HQPGAVFPLG PQEEVLS
Source	E.coli
Target Names	TRIM65
Protein Names	Recommended name: Tripartite motif-containing protein 65
Expression Region	2-517
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.