



Recombinant Mouse GDNF family receptor alpha-like (Gfral), partial

Product Code	CSB-YP750964MO1
Abbreviation	Gfral
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6SJE0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	Q TNDCAHLIQK CLIDANGCEQ SWRSMEDTCL TPGDSCKINN SLHCNLSIQA LVEKNFQFKE CLCMDLHCT VNKLFGKKCT NKTDNMEKDN KDKWNLTTTP FYHGFKQMQS CLEVTEACVG DVVCNAQLAL YLKACSANGN LCDVKHCQAA IRFFYQNMPF NTAQMLAFCD CAQSDIPCQQ SKETLHSKPC ALNIVPPPTC LSVIHTCRND ELCRTHYRTF QTECWPHITG KHEDETCIS MLGKQDLTCS GSESCRAAFL GTFGTVLQVP CACRGVTQAE EHVCMIFQHM LHSKSCFNYP TPNVKDISSY EKKNSKEITL TGFNSFFNG
Source	Yeast
Target Names	Gfral
Protein Names	Recommended name: GDNF family receptor alpha-like
Expression Region	20-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.