



Recombinant Human Ras-like protein family member 11A (RASL11A)

Product Code	CSB-BP750972HU
Abbreviation	RASL11A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6T310
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRPLSMSGHF LLAPIPESSS DYLLPKDIKL AVLGAGRVGK SAMIVRFLTK RFIGDYEPNT GKLYSRLVYV EGDQLSLQIQ DTPGGVQIQD SLPQVVDLSL KCVQWAEGL LVYSITDYDS YLSIRPLYQH IRKVHPDSKA PVIIVGNKGD LLHARQVQTQ DGIQLANELG SLFLEISTSE NYEDVCDVFQ HLCKEVSKMH GLSGERRRAS IIPRPRSPNM QDLKRRFKQA LSPKVKAPSA LG
Source	Baculovirus
Target Names	RASL11A
Protein Names	Recommended name: Ras-like protein family member 11A
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.