



Recombinant Human Zona pellucida-binding protein 2 (ZPBP2)

Product Code	CSB-BP751060HU
Abbreviation	ZPBP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6X784
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	PRFTLFNKKG FIYGKTGQPD KIYVELHQNS PVLICMDFKL SKKEIVDPTY LWIGPNEKTL TGNNRINITE TGQLMVKDFL EPLSGLYTCT LSYKTVKAET QEEKTVKKRY DFMVFAYREP DYSYQMAVRF TTRSCIGRYN DVFFRVLKKI LDSLISDLSC HVEIPSYKCH SVEIPEHGLI HELFIAFQVN PFAPGWKGAC NGSVDCEDTT NHNILQARDR IEDFFRSQAY IFYHNFNRTL PAMHFVDHSL QVVRLDSCR P GFGKNERLHS NCASCCVCS PATFSPDVNV TCQTCVSVLT YGAKSCPQTS NKNQQYED
Source	Baculovirus
Target Names	ZPBP2
Protein Names	Recommended name: Zona pellucida-binding protein 2 Alternative name(s): ZPBP-like protein
Expression Region	21-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.