



# Recombinant Human F-box/WD repeat-containing protein 12 (FBXW12)

<b>Product Code</b>	CSB-BP751062HU
<b>Abbreviation</b>	FBXW12
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6X9E4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEIRLPDLAL KRIFSFLDLF GLLQVSQVNK HWNRIADSDY LWRSLSLQRW DCSNFTNQHL GTHTWKQFFL HQRRKELRLA LAQPHNFIYK VTKNIAFETE LAYLSGNRLT VDEQEKSIC SVSPKQELCA WDVQEGTMIW SSPVQEFHFS NLVTLPQMHL AITMDRKKTI KVVNCQDRDA LAVLPMPQPC YCMEAYLTKD GPFLMVGDAAGDIYTFITLPGLRDVSKVTAFQYGIVLLHCS PDKKWVFCG TYSRTLTPQVF LTESLLRPSE GSVPLSTFLP HKLCASACWT PKVKNRITLM SQSSTGKKTE FITFDLTTKK TGGQTVIQAY EIASFQVAAH LKCPIMWGAS DGYMIVFTSG PYLLLFSITG FLLQRFEDHQ AAINNFVWDP CYVLTTSENS VHVYMWEEGG RHPYLRSCCH LENTWHDHTT DSCISSVMCD NASIVLRVRK VSDSSILVMY SLNT
<b>Source</b>	Baculovirus
<b>Target Names</b>	FBXW12
<b>Protein Names</b>	Recommended name: F-box/WD repeat-containing protein 12 Alternative name(s): F-box and WD-40 domain-containing protein 12 F-box only protein 35
<b>Expression Region</b>	1-464
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.