



Recombinant Human F-box/WD repeat-containing protein 12 (FBXW12)

Product Code	CSB-MP751062HU
Abbreviation	FBXW12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6X9E4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEIRLPDLAL KRIFSFLDLF GLLQVSQV NK HWNRIADSDY LWRSLSLQRW DCSNFTNQHL GTHTWKQFFL HQRRKELRLA LAQPHNFIYK VTKNIAFETE LAYLSGNRLT VDEQEKSII C SVSPKQELCA WDVQEGTMIW SSPVQEFHFS NLVTLPQMHL AITMDRKKTI KVVNCQDRDA LAVLPMPQPC YCMEAYLTKD GPFLMVGDAAGDIYTF TLPGLRDVSKVTAF QYGIVLLHCS PDKKWVFCG TYSRTL PQVF LTESLLRPSE GSVPLSTFLP HKLCASACWT PKVKNRITLM SQSSTGKKTE FITFDLTTKK TGGQTVIQAY EIASFQVA AH LKCP IWMGAS DGYMIVFTSG PYLLLFSITG FLLQRFEDHQ AAINNFVWDP CYVLTTSENS VHVYMWEEGG RHPYLRSCCH LENTWHDHTT DSCISSVMCD NASIVLRVRK VSDSSILVMY SLNT
Source	Mammalian cell
Target Names	FBXW12
Protein Names	Recommended name: F-box/WD repeat-containing protein 12 Alternative name(s): F-box and WD-40 domain-containing protein 12 F-box only protein 35
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.