



Recombinant *Oryza sativa* subsp. *japonica* Germin-like protein 8-2 (GER3)

Product Code	CSB-MP751126OFG
Abbreviation	GER3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YZA9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	SDPSPL QDFCVADLNS PVRVNGFVCK NPMNASADDF FKAAMLDKPR DTNNKVGSNV TLVNVLQLPG LNTLGISIAR LDFAPLGLNP PHTHPRATEI FTVLEGLTYV GFVTSNPDNR LLSKVLNKGD VFVFPEGLIH FQFNPNPHKP AVAIAALSSQ NPGVITIANA VFGSNPPISD DILMKAFQVD KKIIDLQQAQ F
Source	Mammalian cell
Target Names	GER3
Protein Names	Recommended name: Germin-like protein 8-2 Alternative name(s): Germin-like protein 16 Germin-like protein 3 Short name= OsGER3
Expression Region	25-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.