



# Recombinant *Oryza sativa* subsp. *japonica* Hexokinase-4, chloroplastic (HXK4)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP751129OFG   |
| <b>Abbreviation</b>      | HXK4  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q6Z398  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | CSA VAPTSIAPI LADLRLRCAA PLPVLRRVAD AMASGMRAGL<br>ADDGAGELKM IPSHVYSLPT GNETGLFYAL DLGGTNFRVL RVQLGGKDKR<br>IIDTEFEQVS IPREIMHGIT EDLFDFIASG LSRFVATEGD KFHLPPQGRKR<br>ELGFTFSFPV NQTSIDSGIL IKWTKGFAVS GTAGKDVVAC LNAAMERQGL<br>DMRVSALVND TVGTLAGARY WDDVVMVAVI LGTGTNACYI QRTEAIPKLQ<br>HLKLETGNTI INTEWGAFSD GLPLTEFDRE MDDESINPGE QIFEKTISGM<br>YLGEIVRRVL VKMAEVSDLF GHSFPKKLAE PFVLRTPHLC AMQQDTSNLL<br>GEVESILSDV IGVSQASLLA RRVTVEVSDC IIRGGRLAG AGIVGILEKM<br>ENDSRGHIFG RRTVVAMDGG LYEKYPQYRR YMKEAVAELL GPERSNRIAI<br>EHTKDGSGIG AALLAAANSK YAAAQISTR |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | HXK4  |
| <b>Protein Names</b>     | Recommended name: Hexokinase-4, chloroplastic EC= 2.7.1.1 Alternative name(s): Hexokinase II  |
| <b>Expression Region</b> | 38-509  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.