



# Recombinant Human RNA-binding protein MEX3B (MEX3B)

<b>Product Code</b>	CSB-BP751152HU
<b>Abbreviation</b>	MEX3B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6ZN04
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPSSLFADLE RINGSGGGGGG SSGGGETLDD QRALQLALDQ LLLGLDSDE GASLYDSEPR KKSVMTECV PVPSSSEHVAE IVGRQGCKIK ALRAKTNTYI KTPVRGEEPV FVVTGRKEDV AMARREIISA AEHFMSIRAS RNKNTALNGA VPGPPNLPQG TTIQVRVPYR VVGLVVGPKG ATIKRIQQQT HTYIVTPSRD KEPVFEVTGM PENVDREAE IEAHIALRTG GIELTDEND FHANGTDVGF DLHHGSGGSG PGSLWSKPTP SITPTGRKP FSSYRNDSSS SLGSASTDSY FGGGTSSSAA ATQRLADYSP PSPALFAHN GNNNNNGNGY TYTAGGEASV PSPDGCPELQ PTFDPAPAPP PGAPLIWAQF ERSPGGGPAA PVSSSCSSSA SSSASSSSV FPGGGASAPS NANLGLLVHR RLHPGTSCPR LSPPLHMAPG AGEHHLARRV RSDPGGGGLA YAAYANGLGA QLPGLQPSDT SGSSSSSSSS SSSSSSSSGL RRKGSRDCSV CFESEVIAAL VPCGHNLFCM ECANRICEKS EPECPVCHTA VTQAIRIFS
<b>Source</b>	Baculovirus
<b>Target Names</b>	MEX3B
<b>Protein Names</b>	Recommended name: RNA-binding protein MEX3B Alternative name(s): RING finger and KH domain-containing protein 3 RING finger protein 195
<b>Expression Region</b>	1-569
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.