



Recombinant Human RNA-binding protein MEX3B (MEX3B)

Product Code	CSB-BP751152HU
Abbreviation	MEX3B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZN04
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MPSSLFADLE RNGSGGGGGG SSGGGETLDD QRALQLALDQ LSSLGLDSDE GASLYDSEPR KKSVMNTECV PVPSSSEHVAE IVGRQGCKIK ALRAKTNTYI KTPVRGEEPV FVVTGRKEDV AMARREIISA AEHFSMIRAS RNKNTALNGA VPGPPNLPQG TTIQVRVPYR VVGLVVGPKG ATIKRIQQQT HTYIVTPSRD KEPVFEVTGM PENVDREER IEAHIALRTG GIIELTDEND FHANGTDVGF DLHHGSGGSG PGSLWSKPTP SITPTGRKP FSSYRNDSSS SLGSASTDSY FGGGTSSSAA ATQRLADYSP PSPALSFAHN GNNNNNGNGY TYTAGGEASV PSPDGCPELQ PTFDPAPAPP PGAPLIWAQF ERSPGGGPAA PVSSSCSSSA SSSASSSSV FPGGGASAPS NANLGLLVHR RLHPGTSCPR LSPPLHMAPG AGEHHLARRV RSDPGGGGLA YAAYANGLGA QLPGLQPSDT SGSSSSSSSS SSSSSSSSGL RRGSRDCSV CFESEVIAAL VPCGHNLFCM ECANRICEKS EPECVCHTA VTQAIRIFS
Source	Baculovirus
Target Names	MEX3B
Protein Names	Recommended name: RNA-binding protein MEX3B Alternative name(s): RING finger and KH domain-containing protein 3 RING finger protein 195
Expression Region	1-569
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.