



Recombinant *Ashbya gossypii* Enhancer of mRNA-decapping protein 1 (EDC1)

Product Code	CSB-MP751728DOT
Abbreviation	EDC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q753U4
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSSDTMFINS ARLLPDENRK VKQLQRPDRK KHTRGSYAAQ KQQLPNGEQP NFGHGQKQSR KRGSQRQKGR DGAAGDSANV GLTEDLKQLL SIPSGSAGSE SAQKETSAGQ PATGAVADRK RSVPAAGPAG KSSSEPASAS SAVDGMPGLG FQSAHVPFSS PVDNHPMLAE PAATGPVLSV PFPGHVAGYM PCVNNTRVQY PLQTRNNVPM SQPMSQPMSQ PMSQPMSQPM SQPISPYLPQ GFQMHPQSFV GLPVAPMYPQ YMSAQYQSPK QQSACQMPLR AQPPLIPSIT RNALPCVSSA NSTAKASVRS KGSPGSEGGG RRSQNWKSSQ NVGYAGATFA TDQPALSSLP KPSFV
Source	Mammalian cell
Target Names	EDC1
Protein Names	Recommended name: Enhancer of mRNA-decapping protein 1
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.