



Recombinant *Ashbya gossypii* Heterogeneous nuclear rnp K-like protein 2 (HEK2)

Product Code	CSB-YP751746DOT
Abbreviation	HEK2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q756R8
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTMTDSSSPA PEGPAAPAEA AAPAKVDAQA PLVTLRLLLS LKEAARILGR GGNTIDNIRK ANGVKIGIST KEKSCSDRLL EVTGSIDAVA NSLADVVKVL TADDTNEEDA EPEPEQHIFK HLNFILPPPS PDEAEDPMKV KQIGNLRLIV TNSQVSSIIG TAGSKIKKLI DTHSSKLVVS KTFLPDSQDR ILEIQGFNS IANCIKDISQ TLIKDDVLDL KEKRYYPHSK HSKDIHVTAT VAIPAEYVGA LLGHGGNRIA NLRKYTRTKI TVAQEPNQQN EREFTITGND QKSVKLAQNM IMKNLGTEKE RRLEKQ
Source	Yeast
Target Names	HEK2
Protein Names	Recommended name: Heterogeneous nuclear rnp K-like protein 2 Alternative name(s): KH domain-containing protein 1
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.