



Recombinant *Ashbya gossypii* Ribonuclease T2-like (RNY1)

| | |
|--------------------------|---|
| Product Code | CSB-EP751784DOT-B |
| Abbreviation | RNY1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q75BW5 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>) |
| Purity | >85% (SDS-PAGE) |
| Sequence | FQEFITK FPMPMMYNFP SCSSTIPSTC RNETAIADTC CFEYPGGLIL HSQFWNAPYR KRSYRDFGPD DSFTIHGLWN DRCDGSWDQF CRRGSSIRSV VDILSKDSL N RGGLPITGKA LLRQMSMYWK GDRGDENLWV HEY NKHGLCL NTLRPECYQR WGSVASAEDQ AIYDYFRIAM NLHLKIDAYH ALSRQGIKPR CDAPYDAVRM QNALADDFGR EVQMQCTG NR LTGVTTYLL RGGILSENFQ PVDPTQSSSC RGKIYWIPKS GC |
| Source | <i>E.coli</i> |
| Target Names | RNY1 |
| Protein Names | Recommended name: Ribonuclease T2-like Short name= RNase T2-like EC= 3.1.27.1 |
| Expression Region | 24-292 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |