



Recombinant *Ashbya gossypii* ATP-dependent RNA helicase DBP5 (DBP5)

Product Code	CSB-MP751786DOT
Abbreviation	DBP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75C39
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSGNAQKPKDA SEMLADLDLN KTKKETTLEA GSAESAPSAP AENQAKSDSN LINSEYEVKV RLADIQADPN SPLYSVKSFE ELGLAPELLK GLYAMKFQKP SKIQERALPL LLHNPPRNMI AQSQSGTGKT AAFSLTMLSR VDVAVPATQA ICLAPSRELA RQTLVIQEM GKFTKIASQL IVPDSYEKNK AINAHIIVGT PGTVLDMRRL KMIQLGKVKT FVLDEADNML DKQGLGDQCI RVKKFLPKDT QLVLFSAITFD DSVREYARRV VFNANSLELQ RNEVNVSAIK QLFMDCNDER HKFTVLCGLY GLLTIGSSII FVQTKQTANM LYTELKREGH QVSILHGDLQ SADRDRLLIGD FREGRSKVLI TTNVLRGID IPTVSMVVNY DLPMTANGQP DPSTYVHRIG RTGRFRGRTGV AISFIHDKKS YETLAAIQSY FGDIQITKVP TDDMDEMEKI VKKVLK
Source	Mammalian cell
Target Names	DBP5
Protein Names	Recommended name: ATP-dependent RNA helicase DBP5 EC= 3.6.4.13
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.