



Recombinant *Ashbya gossypii* Pre-rRNA-processing protein ESF2 (ESF2)

Product Code	CSB-MP751795DOT
Abbreviation	ESF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75DA3
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSDNDSDVQD FSSEEEDHGL LIDRKKQQTL DFAQDGS DSE GGS DLDDDDI AEEESPIQNE TKAATGEEPA TADGDNIAPQ SADEL RPEER ADQLKQERLR RLK K LKASKK SNHKTGVVYL SKIPPYMKPA KMRQILSRFG DLDRLFLKRE DEHSHRQRVK GGGNKKVMFR EGWAEFIRKK DAKLCAETLN GNIIGGKKGN FYHDDVMNVK YLSGFKWADL TEQIARENDV RQSKLQLEIS QANKLNAEFI RNVEKSKMLN NIRAKKHKDE DSSERTESHL PVKQRKVT SN RASAPESQKQ PASEKLSNVL HNFL
Source	Mammalian cell
Target Names	ESF2
Protein Names	Recommended name: Pre-rRNA-processing protein ESF2 Alternative name(s): 18S rRNA factor 2
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.