



# Recombinant *Ashbya gossypii* AP-3 complex subunit sigma (APS3)

<b>Product Code</b>	CSB-MP751810DOT
<b>Abbreviation</b>	APS3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q75F71
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) ( <i>Eremothecium gossypii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MIHAVLIFNK KCQPRLVKYY TPVDLPKQKL LLEQVYELIS QRNSSIQSSF LITPPSLLSS GSETINEDIQ IYKNYATLY FTFIVDDQES ELAILDLIQT FVEALDRCFA EVNELDLIFN WQTLESVLEE IIQGGMVIET NVKKIVETVD ELNRTSNQEA RFGNGLGNAF QAITMGGFSN WGARQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	APS3
<b>Protein Names</b>	Recommended name: AP-3 complex subunit sigma Alternative name(s): AP-3 complex sigma3A subunit Adapter-related protein complex 3 sigma subunit Clathrin-associated/assembly/adaptor protein, small 3 Sigma-adaptin 3A Sigma3-adaptin
<b>Expression Region</b>	1-185
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.