



# Recombinant *Oryza sativa* subsp. japonica Solanesyl-diphosphate synthase 2, chloroplastic (SPS2)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP751816OFG-B   |
| <b>Abbreviation</b>      | SPS2  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q75HZ9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. japonica (Rice)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | AEAISSAR TTTMIPERIS VSSLLEVSD DLLKLNNNLK SLVGAENPVL<br>VSAAEQIFGA G GKRLRPALV FLVSRATAEL AGLLELTTEH QRLAEIEMI<br>HTASLIHDDV IDDSGMRRGK ETIHQLYGTR VAVLAGDFMF AQSSWFLANL<br>ENIEVIKLIS QVIKDFASGE IKQASTLFDC DVTLLDDYLLK SYYKTASLLA<br>SSTRSAAIFS GVSTTICEQM YEYGRNLGLS FQVVDILDF TQSAEQLGKP<br>AGSDLAKGNL TAPVIFALQD EPKLR EIDS EFSEDSLAT AIDLVHRSGG<br>IRRAQELAKE KGD LALQNLQ CLPKSQFRST LENVVKYNLQ RID    |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | SPS2  |
| <b>Protein Names</b>     | Recommended name: Solanesyl-diphosphate synthase 2, chloroplastic Short name= OsSPS2 EC= 2.5.1.84 Alternative name(s): All-trans-nonaprenyl-diphosphate synthase 2 (geranyl-diphosphate specific)   |
| <b>Expression Region</b> | 63-403  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.