



Recombinant Arabidopsis thaliana Ethylene-responsive transcription factor RAP2-6 (RAP2-6)

Product Code	CSB-EP752043DOA-B
Abbreviation	RAP2-6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7G1L2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MVSMMLTNVVS GETEPSASAT WTMGHKRERE EFSLPPQPLI TGSAVTKECE SSMSLERPKK YRGVRQRPWG KWAAEIRDPH KATRVWLGTG ETAEAAARAY DAAALRFRGS KAKLNFPENV GTQTIQRNSH FLQNSMQPSL TYIDQCPTLL SYSRCMEQQQ PLVGMLQPTE EENHFFEKPW TEYDQYNYSS FG
Source	E.coli
Target Names	RAP2-6
Protein Names	Recommended name: Ethylene-responsive transcription factor RAP2-6 Alternative name(s): Protein RELATED TO APETALA2 6
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.