



Recombinant Drosophila melanogaster Protein N-terminal glutamine amidohydrolase (tun)

Product Code	CSB-BP752082DLU
Abbreviation	tun
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7K2Y9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MTTDFLFPKI SDCSYVSCYC EENVWKLCEQ VKRTRPEELS KCYAVFVSNE GRTVPLWRQK AGRGDDQVVI WDYHVFFIHN PLLNRCLVFD LDTTLFPPTY FHKYVTETFR SDLALRPEHH RFFRVIPADT YLIEFSSDRR HMRRPDGSWI KPPPSYPPIL SNSNMHCLGD FICMSAGKGP GSVYSLSEFV QNFYKSPHVM AQNNK
Source	Baculovirus
Target Names	tun
Protein Names	Recommended name: Protein N-terminal glutamine amidohydrolase EC= 3.5.1.- Alternative name(s): Protein NH2-terminal glutamine deamidase Short name= N-terminal Gln amidase Short name= Nt(Q)-amidase Protein tungus
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.