



Recombinant Human Basic leucine zipper and W2 domain-containing protein 1 (BZW1)

Product Code	CSB-MP752097HU
Abbreviation	BZW1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L1Q6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MNNQKQQKPT LSGQRFKTRK RDEKERFDPT QFQDCIIQGL TETGTDLEAV AKFLDASGAK LDYRRYAETL FDILVAGGML APGGTLADDM MRTDVCVFAA QEDLETMQAF AQVFNKLIRR YKYLEKGFED EVKLLLLFLK GFSESERNKL AMLTGVLLAN GTLNASILNS LYNENLVKEG VSAAFVAVKLF KSWINEKDIN AVAASLRKVS MDNRLMELFP ANKQSVEHFT KYFTEAGLKE LSEYVRNQQT IGARKELQKE LQEQMSRGDP FKDIILYVKE EMKKNNIPEP VVIGIVWSSV MSTVEWNNKKE ELVAEQAIKH LKQYSPLLA FTTQGQSELT LLLKIQEYCY DNIHFMKAFQ KIVVLFYKAE VLSEEPILKW YKDAHVAKGK SVFLEQMKKF VEWLKNAEEEE SESEAEEGD
Source	Mammalian cell
Target Names	BZW1
Protein Names	Recommended name: Basic leucine zipper and W2 domain-containing protein 1 Alternative name(s): Protein Orf
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.