



Recombinant Human ATP-binding domain-containing protein 4 (ATPBD4)

Product Code	CSB-EP752102HU-B
Abbreviation	ATPBD4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7L8W6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRVAALISGG KDSCYNMMQC IAAGHQIVAL ANLRPAENQV GSDELDSYMY QTVGHHAILD YAEAMALPLY RRTIRGRSLD TRQVYTKCEG DEVEDLYELL KLVKEKEEVE GISVGAILSD YQRIRVENVC KRLNLQPLAY LWQRNQEDLL REMISSNIQA MIIKVAALGL DPKHLGKTL DQMEPYLIEL SKKYGVHVCG EGGEYETFTL DCPLFKKKII VDSSEVVIHS ADAFAPVAYL RFLELHLEDK VSSVPDNYRT SNYIYNF
Source	E.coli
Target Names	DPH6
Protein Names	Recommended name: ATP-binding domain-containing protein 4
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.