



Recombinant Danio rerio Protein FAM172A (fam172a)

Product Code	CSB-MP752606DIL
Abbreviation	fam172a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7T297
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	QL LAQDEKEKTT ALKDLLSRID LDELMKKDEP PLDFPETLDG FEYAFNEQQG LRHTKTGEPF VFNYKEDLHR WNQKRYEALG EIITKYVYER LEKDCNLKKE ILPVDATEDE PKSFIYVSED ALSNPDKIMV LIQGGSGVVRA GQWARRLIIN EDLNSGTQIP FINRAKEEGY GVIVLNPEN SLEVEKVSDP AKKASPDASD EPAEKRERKE ERENKKKKDF YEKYRNPQKE REMEHIPIRE NGSPEEHTLY AWDYFISKSL AKNIFIVAHS YGGLSFVNLN INREAEVKSR VRAVAMTDSV HNVWHQEASK SIQDWLREHC CNWVSSPEPL DTLVDPMLND CPRVSAGTEK HELTSWKSFQ SIFKFFSNAL KDHEDDSSSS PVTTRSSHRI KHEDF
Source	Mammalian cell
Target Names	fam172a
Protein Names	Recommended name: Protein FAM172A
Expression Region	19-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.