



Recombinant Mouse Wilms tumor protein 1-interacting protein (Wtip)

Product Code	CSB-EP752651MO
Abbreviation	Wtip
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7TQJ8
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MQRSRTAADD AALLLAGLGL RESEPTAGSP GRVRRGPRAV DEAAPASGRR GKGGCGGPEA APDVPSRPER GPRASLAGSD GGSARSSGIS LGYDQRHGPG PGPPSGGSAR SSVSSLGSRG SAGACADLLP PGVGPAPARS PEPAQFPFPL PSLPLPPGRE GGPSAAERRL EALTRELER LEARTARDYF GICIKCGLGI YGARQACQAM GSLYHTDCFI CDSCGRRLRG KAFYNVGEKV YCQEDFLYSG FQQTADKCSV CGHLIMEMIL QALGKSYHPG CFRCSVCNEC LDGVPFTVDV DNNIYCVRDY HTVFAPKCAS CARPILPAQG CETTIRVSM DRDYHVECYH CEDCGLQLSG EEGRRCYPLE GHLLCRRCHL RRLGQGPLPS PAVHVTEL
Source	E.coli
Target Names	Wtip
Protein Names	Recommended name: Wilms tumor protein 1-interacting protein Short name= WT1-interacting protein
Expression Region	1-398
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.