



# Recombinant Human Actin-binding LIM protein 2 (ABLIM2)

<b>Product Code</b>	CSB-BP753612HU
<b>Abbreviation</b>	ABLIM2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6H8Q1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAVSQPQAA PSPLEKSPST AILCNTCGNV CKGEVLRVQD KYFHIKCFVC KACGCDLAEG GFFVRQGEYI CTLDYQRLYG TRCFSCDQFI EGEVVSALGK TYHPDCFVCA VCRLPFPPGD RVTFNGKECM CQKCSLPVSV GSSAHLSQL RSCGGCGTEI KNGQALVALD KHWHLGCFKC KSCGKLLNAE YISKDGLPYC EADYHAKFGI RCDSCEKYIT GRVLEAGEKH YHPSCALCVR CGQMFAEGEE MYLQGSSIIWH PACRQAARTE DRNKETRTSS ESIISVPASS TSGSPSRVIY AKLGGEILDY RDLAALPKSK AIYDIDRPDM ISYSPYISHS AGDRQSYGEG DQDDRSYKQC RTSSPSSTGS VSLGRYTPTS RSPQHYSRPG SESGRSTPSL SVLSDSKPPP STYQQAPRHF HVPDTGVKDN IYRKPPIYRQ HAARRSDGED GSLDQDNRKK SSWLMLKGDA DTRTNSPDLT TQSLSHSSGT DRDPLQRMAG DSFHSRFPYS KSDPLPGHGK NGLDQRNANL APCGADPDAS WGMREYKIYP YDSLIVTNRI RVKLPKDVD RTRLERHLSPE EFQEVFGMSI EEFDRALWK RNDLKKKALL F
<b>Source</b>	Baculovirus
<b>Target Names</b>	ABLIM2
<b>Protein Names</b>	Recommended name: Actin-binding LIM protein 2 Short name= abLIM-2 Alternative name(s): Actin-binding LIM protein family member 2
<b>Expression Region</b>	1-611
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.