



# Recombinant Rat Developmental pluripotency-associated protein 3 (Dppa3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP753700RA-B  |
| <b>Abbreviation</b>      | Dppa3   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q6IMK0  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Rattus norvegicus (Rat)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDEPSEKVDV VVNPETQMYD GSQREDEGDS PDDSEILQPE TLVKVMKKLT<br>LNPSAKPTKY HRRQRVRLQV KSQPVENRSE RIMREVQSAF PRRRVRTLLS<br>VLKDPIARMR RFVRIEQRQR QLEGNERRDE PFRCLCTFCH YQRWDPSENA<br>KIGQNQKN  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Dppa3   |
| <b>Protein Names</b>     | Recommended name: Developmental pluripotency-associated protein 3   |
| <b>Expression Region</b> | 1-158   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This gene encodes a protein that in mice may function as a maternal factor during the preimplantation stage of development. In mice, this gene may play a role in transcriptional repression, cell division, and maintenance of cell pluripotentiality. In humans, related intronless loci are located on chromosomes 14 and X.   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |