



Recombinant Human Protein FAM166A (FAM166A)

Product Code	CSB-BP753727HU
Abbreviation	FAM166A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6J272
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTTTQKHDLF TPEPHYVPGY AGFFPQLRYQ VGNTYGRRTTG QLLTDPSVQK SPCSVLSPMS KPKFIEDFSQ SKPPRVPCQD LTEPYIPHYT SLKPSKNFEI LGQLPPLEVD AQEPPGVENI PRQILLPAGF TPDTPHPPCP PGRKGDSRDL GHPVYGEEAW KSATPVCEAP RQHQLYHCQR DEYPPPARRQ QETLDVGSFQ RLPQLDHPNL IQRKAISGYA GFIPRFTWVM GLNYRDGVMQ AMDEFDKSQF LFRNPHCDLG EKLPGTHWPS NHIYSSQGLI PFYMGFIPAM QDNYALTFGN STRRAYWKEW AKRNHTL
Source	Baculovirus
Target Names	FAM166A
Protein Names	Recommended name: Protein FAM166A
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.