



Recombinant *Oryza sativa* subsp. japonica

Putative cyclin-F2-1 (CycF2-1)

Product Code	CSB-YP753743OFG
Abbreviation	CycF2-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K1Z6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSIMEPYVA DLLADDITAS MVELLSGDGG AAQMDVGVLD AYLRAIGALP AHPAAPGADL AAAAEVESMA SNDDTNGNWD TKVDAKVP SA FLPPPPGFPP LPVPALADEP VYAAPVDEGD AIRAFMQQLE WSEQYNGDDD APAPDDSMAS RPQLCAPYDD DIDANLRAME KDAAERPSPD YLDTVHNGQI SAASRASLVA WMGRLTHRYE LAAGTLHRAV SYFDRFLSAR ALPSYTEHQL SLVGATAVYT AAKYEDQGT VFKLDAREIAS YGEFASAQEV LAMEREMMAA LGYRLGGPNA ETFVEHFTRY SKGKEELRVQ RLARHIADRS LESYGCLGYL PSVVAAAVIS IARWTLNPPG ALPWSSELHE LTGYSSQDIS SCVLTVLNTQ
Source	Yeast
Target Names	CycF2-1
Protein Names	Recommended name: Putative cyclin-F2-1 Short name= CycF2;1
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.