



Recombinant *Oryza sativa* subsp. japonica Probable plastid-lipid-associated protein 2, chloroplastic (PAP2)

Product Code	CSB-MP753744OFG
Abbreviation	PAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K439
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	AVAGDA EDEWGKEPAA DQGGAAAAVA EAPADVPTS EVAELKAKLK EALYGTERGL RASSETRAEV VELITQLEAR NPTPAPTEAL TLLNGKWILA YTSFSQLFPL LGSGSLPQLV KVEEISQTID SENFTVQNCI KFSGPLATTS VSTNAKFEVR SPKRQVQIKFD EGIIGTPQLT DSIVLPEKFE LFGQNIDLTP LKGIFSSIEN AASSVARTIS GQPPLKIPIR TDNAESWLLT TYLDDELRLS RGDGSSIFVL FKEGSTLLY
Source	Mammalian cell
Target Names	PAP2
Protein Names	Recommended name: Probable plastid-lipid-associated protein 2, chloroplastic Alternative name(s): Fibrillin-like protein 2
Expression Region	55-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.