



Recombinant *Oryza sativa* subsp. japonica Putative germin-like protein 2-3 (Os02g0491800, LOC_Os02g29020)

Product Code	CSB-EP753745OFG-B
Abbreviation	Os02g0491800, LOC_Os02g29020
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6K5P8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	SD PSHLQDLCVA DKASTVRVNG VACKDGEDVA AEDFFFSGLH MAGNTTNKQG SAVTAVNVAQ VPGLNTLGIS LARIDYALHG LNPPHHPRA TEILTVLEGS LYVGFVTSNP ENKLFTKVIN KGDVVFVFPKG LVHFQFNRYGT TDAVAIVAL S QNPGVITVA NAVFGSKPSI TDDILAKAFQ VEKTVVDQIQ AKF
Source	E.coli
Target Names	Os02g0491800
Protein Names	Recommended name: Putative germin-like protein 2-3
Expression Region	29-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.