



# Recombinant *Oryza sativa* subsp. japonica Ubiquitin-like modifier-activating enzyme 5 (Os02g0506500, LOC\_Os02g30310)

<b>Product Code</b>	CSB-EP753746OFG
<b>Abbreviation</b>	Os02g0506500, LOC_Os02g30310
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6K6K7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDEEQLRALL RDL DALKQRP DPAAIDRMRE RVAGMVTPAA AARSKIKDMS SEVVDSPNYS RLMALQRMGI VDNYERIRDY SIAIVGIGGV GSVAEMLTR CGIGRLLLYD YDTVELANMN RLFFRPDQVG MTKTDAAVQT LSGINPDVTL ESYSLNITTV KGFETFLGSL KARSSDGRNT GVDLVLSCVD NYEARMVVNQ ACNELGQWTM ESGVSEDAVS GHIQLLVPGE TACFACAPPL VVASGVDERT LKREGVCAAS LPTTMGVVAG LLVQNALKYL LKFGQVSPYL GYNLKDYFP TMEKPNPQC SNPACVQRQK EYMQSKPARD AAKAKMEAE ASAADECPVH LDNDWNISVV DDSDTVTPSI LSTGADSLPE GLVRELPTAD SYQEPVAPVT SGAIDDDLEE LQRQLDALNS S
<b>Source</b>	E.coli
<b>Target Names</b>	Os02g0506500
<b>Protein Names</b>	Recommended name: Ubiquitin-like modifier-activating enzyme 5 Short name= Ubiquitin-activating enzyme 5
<b>Expression Region</b>	1-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.