



Recombinant Human HLA class I histocompatibility antigen protein P5 (HCP5)

Product Code	CSB-EP754020HU
Abbreviation	HCP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6MZN7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLLRMSEHRN EALGNYLEMR LKSSFLRGLG SWKSNPLRLG GWTILLTLTM GQGEPGGPQG DPWVPHELLL PSLCDSSHAS SWGSGSITCA WRGGDSSSHP LVSGHILSNS PVAAVMCSSM GTHLSPFKGT LL
Source	E.coli
Target Names	HCP5
Protein Names	Recommended name: HLA class I histocompatibility antigen protein P5 Alternative name(s): HLA complex protein P5 Protein P5-1
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.