



Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 65 (At5g01700)

Product Code	CSB-YP754100DOA
Abbreviation	At5g01700
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NKS1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGVCCSKGTG IIVEHGADDG NECGDGAEV RDTNDGAVVR TRGSSKHVSM SIKQGKKGIN QDAMTVWENF GGEEDTIFCG VFDGHGPMGH KISRHCENL PSRVHSKIRS SKSAGDENIE NNSSQSQEEL FREFEDILVT FFKQIDSELG LDSPYDSFCS GTTAVTVFKQ ADCLVIANLG HSRAVLGTRS KNSFKAVQLT VDLKPCVQRE AERIVSCKGR VFAMEEEDPDV YRVWMPDDDC PGLAMSRAGF DFCLKDYGLV CIPDVFCRKV SREDEFVULA TDGIWDVLSN EEVVKVVGSC KDRSVAEML VQRAARTWRT KFPASKADDC AVVVLYLNHR PYPREGNVS R AISTISWRSN KSNNECYGAA PLSPLGLSQR VS
Source	Yeast
Target Names	At5g01700
Protein Names	Recommended name: Probable protein phosphatase 2C 65 Short name= AtPP2C65 EC= 3.1.3.16
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.