



Recombinant Arabidopsis thaliana F-box protein At2g26850 (At2g26850)

Product Code	CSB-EP754101DOA
Abbreviation	At2g26850
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NLB1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLLYLLITCL SFFFFSKSLS FPQWASKTKN LLSFYFSKYL FTNSLHQITP DLASPVLGKM SILDLPDLPL DCILELLPPS ELCTMARVCS SLRERCVSDH LWEKHLKTKW GKILGPSAHK EWQCYLSSPY HLDSPHHKTS HLGLAKIISL MRSLSIFRD DDHRRRYPS IPLDSTMNFY LSLETGRFWF PAQVYNRENG HVGFMSCYD AELSYDTDN TFQARYPPHG KRAIAVEKDV TWERIRAAPV DASPHNLYVS DSLNELKPGD HIEIQWRRNK EFPYGWWYSV VGHMESC DGN LNHCQCHNSE MMVLEFNQYT VGSRWRK TMI NRRDHREKGN EEDGFYGGIR KINCKEDIEM WKRLWPSSIL E
Source	E.coli
Target Names	At2g26850
Protein Names	Recommended name: F-box protein At2g26850
Expression Region	1-371
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.