



Recombinant Human Homeobox-containing protein 1 (HMBOX1)

| | |
|--------------------------|--|
| Product Code | CSB-MP754126HU |
| Abbreviation | HMBOX1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6NT76 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLSSFPVVLL ETMSHYTDEP RFTIEQIDLL QRLRRTGMTK HEILHALETL DRLDQEHSDK FGRRSSYGGS SYGNSTNNVP ASSSTATAST QTQHSGMSPS PSNSYDTSPQ PCTTNQNGRE NNERLSTSNQ KMSPTRYHAN SMGQRSYSFE ASEEDLDVDD KVEELMRRDS SVIKEEIKAF LANRRISQAV VAQVTGISQS RISHWLLQQG SDLSEQKKRA FYRWYQLEKT NPGATLSMRP APIPIEDPEW RQTPPPVSAT SGTFRLLRGS RFTWRKECLA VMESYFNENQ YPDEAKREEI ANACNAVIQK PGKKLSDLER VTSLKVYNWF ANRRKEIKRR ANIEAAILES HGIDVQSPGG HSNSDDVDGN DYSEQDDSTS HSDHQDPISL AVEMA AVNHT ILALARQGAN EIKTEALDDD |
| Source | Mammalian cell |
| Target Names | HMBOX1 |
| Protein Names | Recommended name: Homeobox-containing protein 1 |
| Expression Region | 1-420 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |