



Recombinant *Xenopus tropicalis* Charged multivesicular body protein 2b (chmp2b)

Product Code	CSB-EP754144XBF
Abbreviation	chmp2b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NVL7
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MASLFKKKTV DDIREQNKE LRGTQRAITR DRAALEKQEK QLEMEIKKMA KAGNKDACRV LAKQLVQLRK QKTRTYAVSS KVTSMSTQTK VMSSQMKMAG AMSTTAKTMQ AVNKKMDPQK TLQTMQNFQK ENMKMEMTEE MINDTLDDIF DASDDEEESQ DIVNQVLDEI GIEISGKMAK APSAAKGLPS TSASKSSGIS DEEIERQLKA LGVD
Source	E.coli
Target Names	chmp2b
Protein Names	Recommended name: Charged multivesicular body protein 2b Alternative name(s): Chromatin-modifying protein 2b Short name= CHMP2b
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.