



Recombinant Danio rerio UPF0696 protein C11orf68 homolog (P5436, zgc:63470)

Product Code	CSB-BP754149DIL
Abbreviation	P5436, zgc:63470
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6NWH0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAAESPASP SELSECFSAE DYAAESLAAD LDPWITFDAR KTPRAEFSSW LQSHRPSSVS RYGDGEAGAG PVGWI AVRGP DYSEQSGDVE GLQDSWETLL TSGRSVSFQN IRELALNHSV LDGKWLMLHD TGFKVDRAWE SIARATVLEG KIWSAKVSPR DPSSNSRHVI CVYNQNFTDE EQVVRLDSAI RAAGVKCVLY YKPDVYTYLG IYRNNRWKLC PTIYESMFDL ESVPRRSHIL NKVTNLELS
Source	Baculovirus
Target Names	P5436
Protein Names	Recommended name: UPF0696 protein C11orf68 homolog
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.