



Recombinant Human Protein FAM117B (FAM117B)

Product Code	CSB-EP754170HU-B
Abbreviation	FAM117B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P1L5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSQRVRRNGS PTPAGSLGGG AVATAGGPGS RLQPMRATVP FQLKQQQQQQ HGSPTRRGGG GGGNNGGCC GGASGPAGGG GGGGPRTASR STSPTRGGGN AAARTSPTVA TQTGASATST RGTSPTRSAA PGARGSPPRP PPPPLLGTV SSPSSSPHL WTGEVSAAPP PARVRHRRRS PEQSRSSPEK RSPSAPVCKA GDKTRQPSSS PSSIIIRTSS LDTLAAPYLA GHWPRDSHGQ AAPCMRDKAT QTESAWAEEY SEKKKGSHKR SASWGSTDQL KEIAKLRQQL QRSKHSSRHH RDKERQSPFH GNHAAINQCQ APVPKSALIP VIPITKSTGS RFRNSVEGLN QEIEIIKET GEKEEQLIPQ DIPDGHRAAPP PLVQRSSSTR SIDTQTPGGA DRGSNNSSRS QSVSPTSFLT ISNEGSEESP CSADDLLVDP RDKENGNNSP LPKYATSPKP NNSYMFKREP PEGCERVKVF EECSPKQLHE IPAFYCPDKN KVNFIKSGS AFCLVSILKP LLPTPDLTK GSGHSLTVTT GMTTLLQPI AVASLSTNTE QDRVSRGTST VMPSASLLPP PEPIEEAEG
Source	E.coli
Target Names	FAM117B
Protein Names	Recommended name: Protein FAM117B Alternative name(s): Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein
Expression Region	1-589
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.