



Recombinant Human WD repeat-containing protein 93 (WDR93)

Product Code	CSB-EP754176HU-B
Abbreviation	WDR93
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P2C0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSFPRGSQTQ KIKHPIGTRK GPLEVPPPTE KDWPKDDEQD HVLVDPDEEL DSLQPQPYRMI NKLVNLLFDQ SWEIIEERNA LREAESSQIQ PTVYPPPLGEI QLNKMPNCMA VSQDYVFIGG AKGFSIYNLY SAKQIYAWEK LKVDVTSIWA TDLGNEILIA PVDEMGIIRL FYFYKEGLYL VKAINEVDDT SKQTTCIKME ISQGGDFAAF LLQGAGDIWL DVYKLPKETW LKKLEHPQLT SNPKKKVRQP QLNSLGPISA DPLEMDANVS FKGDIKLSLP VYIMKIKPPK PVTGTTFKSP LEVFAKIKDC YGLGSGQNHF IKDSQWEQQA EIFNASYKKY LDREWEEEP STATFYFLLP SCLFAMPPEV KGPSGMACVL GIHWTRSHNF FLYSLNRTLK DKADPEGVWP CAPIAVSQL SCSSSYLVLA CEDGVLTLWD LAKGFPLGVA ALPQGCFCQS IHFLKYFSVH KGQNMYPEGQ VKSQMKCVVL CTDASLHLVE ASGTQGPTIS VLVERPVKHL DKTICAVAPV PALPGMVLIF SKNGSVCLMD VAKREIICAF APPGAFPLEV PWKPVFAVSP DHPCFLLRGD YSHETASTDD AGIQYSVFYF NFEACPLEN ISKNCTIPQR DLDNMAFPQA LPLEKRCERF LQKSYRKLEK NPEKEEEHWA RLQRYSLSLQ RENFKK
Source	E.coli
Target Names	WDR93
Protein Names	Recommended name: WD repeat-containing protein 93
Expression Region	1-686
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.