



# Recombinant Human WD repeat-containing protein 93 (WDR93)

<b>Product Code</b>	CSB-MP754176HU
<b>Abbreviation</b>	WDR93
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q6P2C0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSFPRGSQTQ KIKHPIGTRK GPLEVPPPTE KDWPKDDEQD HVLVDPDEEL DSLQPQPYRMI NKLVNLLFDQ SWEIIEERNA LREAESSQIQ PTVYPPPLGEI QLNKMPNCMA VSQDYVFIGG AKGFSIYNLY SAKQIYAWEK LKVDVTSIWA TDLGNEILIA PVDEMGIIRL FYFYKEGLYL VKAINEVDDT SKQTTCIKME ISQGGDFAAF LLQGAGDIWL DVYKLPKETW LKKLEHPQLT SNPKKKVRQP QLNSLGPISA DPLEMDANVS FKGDIKLSLP VYIMKIKPPK PVTGTTFKSP LEVFAKIKDC YGLGSGQNHF IKDSQWEQQA EIFNASYKKY LDREWEEEL STATFYFLLP SCLFAMPPEV KGPSGMACVL GIHWTRSHNF FLYSLNRTLK DKADPEGVWP CAPIAVSQL SCSSSYLVLA CEDGVLTLWD LAKGFPLGVA ALPQGCFCQS IHFLKYFSVH KGQNMYPEGQ VKSQMKCVVL CTDASLHLVE ASGTQGPTIS VLVERPVKHL DKTICAVAPV PALPGMVLIF SKNGSVCLMD VAKREIICAF APPGAFPLEV PWKPVFAVSP DHPCFLLRGD YSHETASTDD AGIQYSVFYF NFEACPLLEN ISKNCTIPQR DLDNMAFPQA LPLEKRCERF LQKSYRKLEK NPEKEEEHWA RLQRYSLSLQ RENFKK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WDR93
<b>Protein Names</b>	Recommended name: WD repeat-containing protein 93
<b>Expression Region</b>	1-686
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.