



Recombinant *Xenopus laevis* Serine/threonine-protein kinase 4 (stk4)

Product Code	CSB-BP754221XBE
Abbreviation	stk4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PA14
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	METVQLRNNP RRHLKKLSEE SLNKQPEEVF DVLEKLGEGS YGSVYKASHK ETSQIVAIKQ IPVESDLQEI IKEISIMQQC DSPHVVKYYG SYFKNTDLWI VMEFCGGGSV SDIIRLRKQT LNEDEIATIL QSTLKGLEYL HFMRKIHRDI KAGNILLSCE GTAKLADFGV AGQLTDTMAK RNTVIGTPFW MAPEVIQEI YNCVADIWSL GITAIEMAEG KPPYAEIHPM RAIFMIPSNP PPTFRKPELW SKDFVDFINL CLVKNPELRS SATELLQHPF IKAAGGESIL RHLLNAAQDE KLKRTELKQR EVEPEEKENV NEDEVVGTM VQAGGKDLNT MKELGMMSEG ADGTMVEKDK LETQMGTMLI NDEDEEEETG TMKQCTEPAQ QAKPSFLEYF EQKENQFGTP EKTSPSTDP SEWKIPLNGD YSFLKDWSVA EVQLKLNSLD PMMEREIEEI HHKYQAKRQP ILEAIESKKR RQQNF
Source	Baculovirus
Target Names	stk4
Protein Names	Recommended name: Serine/threonine-protein kinase 4 EC= 2.7.11.1 Cleaved into the following 2 chains: 1. Serine/threonine-protein kinase 4 37kDa subunit Short name= 2. MST1/N 3. Serine/threonine-protein kinase 4 18kDa subunit
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.