



Recombinant Danio rerio Putative fibroblast growth factor 1 (fgf1)

Product Code	CSB-MP754234DIL
Abbreviation	fgf1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PBT8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTEADIAVKS SPRDYKKLTR LYCMNGGFHL QILADGTVAG AADENTYSIL RIKATSPGVV VIEGSETGLY LSMNEHGKLY ASSLVTDSEY FLEKMEENHY NTYQSQKYGE NWYVGIKKNG KMKRGPRTI GQKAIFLPR QVEQEED
Source	Mammalian cell
Target Names	fgf1
Protein Names	Recommended name: Putative fibroblast growth factor 1 Short name= FGF-1 Alternative name(s): Acidic fibroblast growth factor Short name= aFGF Heparin-binding growth factor 1 Short name= HBGF-1
Expression Region	1-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.