



Recombinant Danio rerio 40S ribosomal protein S3a (rps3a)

Product Code	CSB-EP754235DIL
Abbreviation	rps3a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PBY1
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	AVGKNKRLT KGGKKGAKKK IVDPFSKKDW YDVKAPAMFN IRNLGKTLVT RTQGTRIASD GLKGRVFEVS LADLQNDEVA FRKFKLVTEQ VQGKNCLTNF HGMDLTRDKM CSMVKKWQTM IEAHVDVKTG DGYLLRLFCV GFTKKRTNQT RKTSYAQHQQ VRQIRKKMME IMTREVQTND LKEVVNKLIP DSVGKDIEKA CQSIYPLHDV YVRKVKMLKK PKFELGKLME LHGEGGSTSA AAKPAEGDTG AKVERADGYE PPIQESV
Source	E.coli
Target Names	rps3a
Protein Names	Recommended name: 40S ribosomal protein S3a
Expression Region	2-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.