



Recombinant Human F-box/LRR-repeat protein 19 (FBXL19)

Product Code	CSB-MP754242HU
Abbreviation	FBXL19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PCT2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGMKVPKGGE SGPSALLTPP MSSSSRGPGA GARRRRTRCR RRCRACVRTEC GDCHFCDMK KFGGPGRMKQ SCLLRQCTAP VLPHTAVCLL CGEAGKEDTV EGEEEFGLS LMECTICNEI VHPGCLKMGK AEGVINAIEIP NCWECPRCTQ EGRTSKDSGE GPGRRRADNG EEGASLGSGW KLTEEPPLPP PPPRRKGPLP AGPPPEDVPG PPKRKEREAG NEPPTPRKKV KGGRRERHLKK VGGDACLLRG SDPGGPGLLP PRVLNPSQAF SSCHPGLPPE NWEKPKPLA SAEGPAVPSP SPQREKLERF KRMCQLLERV PDTSSSSSDS DSDSDSSGTS LSEDEAPGEA RNGRRPARGS SGEKENRGGR RAVRPGSGGP LLSWPLGPAP PPRPPQLERH VVRPPRSPE PDTLPLAAGS DHPLPRAAWL RVFQHLGPRE LCICMRVCRT WSRWCYDKRL WPRMDLSRRK SLTPPMLSGV VRRQPRALDL SWGVSQKQL MWLLNRLQGL QELVLSGCSW LSVSALGSAP LPALRLLDLR WIEDVKDSQL RELLPPPDPT KPGQTESRGR LQGVAELRLA GLELTDASLR LLLRHAPQLS ALDLSHCAHV GDPSVHLLTA PTSPLRETLV HLNLAGCHRL TDHCLPLFRR CPRLRRLDLR SCRQLSPEAC ARLAAAGPPG PFCRCPEEKLL LKDS
Source	Mammalian cell
Target Names	FBXL19
Protein Names	Recommended name: F-box/LRR-repeat protein 19 Alternative name(s): F-box and leucine-rich repeat protein 19
Expression Region	1-694
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.