



Recombinant Human Kelch-like protein 35 (KLHL35)

Product Code	CSB-MP754256HU
Abbreviation	KLHL35
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6PF15
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRWVRHDAPA RRGQLRRLE HVRLPLLAPA YFLEKVEADE LLQACGECRP LLEARACFI LGREAGALRT RPRRFMDLAE VIVVIGGCDR KGLLKLPFAD AYHPESQRWT PLPSLPGYTR SEFAACALRN DVYVSGGHIN SHDVWMFSSH LHTWIKVASL HKGRWRHKMA VVQGQLFAVG GFDGLRRLHS VERYDPFSNT WAAAAPLPEA VSSAAVASCA GKLFVIGGAR QGGVNTDKVQ CFDPKEDRWS LRSPAPFSQR CLEAVSLEDT IYVMGGLMSK IFTYDPGTDV WGEAAVLSP VESCGVTVCD GKVHILGGRD DRGESTDKVF TFDPSSGQVE VQPSLQRCTS SHGCVTIIQS LGR
Source	Mammalian cell
Target Names	KLHL35
Protein Names	Recommended name: Kelch-like protein 35
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.