



Recombinant Outer membrane protein Omp38 (omp38)

Product Code	CSB-YP754340ABS
Abbreviation	omp38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6RYW5
Product Type	Recombinant Protein
Immunogen Species	Acinetobacter baumannii
Purity	>85% (SDS-PAGE)
Sequence	A NAGVTVTPLL LGYTFQDSQH NNGGKDG NLT NGPELQDDLF VGAALGIELT PWLGF AEY N QVKGDVDGAS AGAEYKQKQI NGNFYVTS DL ITKNYDSKIK PYVLLGAGHY KYDFDGVNRG TRGTSEEGTL GNAGVGAFWR LNDALSLRTE ARATYNADEE FWNYTALAGL NVVLGGHLKP AAPVVEVAPV EPTPVAPQPQ ELTEDLN MEL RVFFDTN KSN IKDQYKPEIA KVAEKLSEYP NATARIEGHT DNTGPRKLNE RLSLARANSV KSALVNEYNV DASRLSTQGF AWDQPIADNK TKEGRAMNRR VFATITGSRT VVVQPGQEAA APAAAQ
Source	Yeast
Target Names	omp38
Protein Names	Recommended name: Outer membrane protein Omp38 Alternative name(s): Outer membrane protein OmpA Outer membrane protein OmpAb
Expression Region	20-356
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.