



Recombinant Human HHIP-like protein 2 (HHIPL2)

Product Code	CSB-BP754409HU
Abbreviation	HHIPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UWX4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	HPQCLDYGP PFQPPLHLEF CSDYESFGCC DQHKDRRIAA RYWDIMEYFD LKRHELCDY IKDILCQEC PYAAHLYDAE NTQTPLRNLP GLCSDYCSAF HSNCHSAISL LTNDRGLQES HGRDGTRFCH LLDLPDKDYC FPNVLRNDYL NRHLGMVAQD PQGCLQLCLS EVANGLRNPV SMVHAGDGTH RFFVAEQVGV VVVYLPDGSR LEQPFLDLKN IVLTPWIGD ERGFLGLAFH PKFRHNRKfy IYYSCLDKKK VEKIRISEMK VSRADPNKAD LKSERVILEI EEPASNHNGG QLLFGLDGYM YIFTGDGGQA GDPFGLFGNA QNKSSLLGKV LRIDVNRAGS HGKRYRVPSP NPFVSEPGAH PAIYAYGIRN MWRCVDRGD PITRQGRGRI FCGDVGQNR EVDLILKGG NYGWRAKEGF ACYDKKLCHN ASLDDVLPYI AYGHAVGKSV TGGYVYRGCE SPNLNGLYIF GDFMSGRLMA LQEDRKNKKW KKQDLCLGST TSCAFPGLIS THSKFIISFA EDEAGELYFL ATSYPsAYAP RGSYKFVDP SRRAPPgKCK YKVPVVRTKS KRIPFRPLAK TVLDLLKEQS EKAARKSSSA TLASGPAQGL SEKGSSKLA SPTSSKNTLR GPGTKKKARV GPHVRQGKRR KSLKSHSGRM RPSAEQKRAG RSLP
Source	Baculovirus
Target Names	HHIPL2
Protein Names	Recommended name: HHIP-like protein 2
Expression Region	42-724
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.