



Recombinant *Oryza sativa* subsp. japonica

Uncharacterized protein

Os08g0218700/LOC_Os08g12160

(Os08g0218700, LOC_Os08g12160)

Product Code	CSB-MP754504OFG
Abbreviation	Os08g0218700, LOC_Os08g12160
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6YUE5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	AKEEKPVVVL PPAAAPGEAP SADAAAFVRS CCDTALQADR DGSSFCYYHL LPYAAFFEGN QVKVAEVAAT ILSTNLWVYV DQLRKVQGGA GKGDPNLNAC VDDFSVAAGE NITREALQSL GRLAAAGNGK RSKEDLENAQ KWIKGVEKPY NGGIGKASGC EIGYLFYSD DLPAQKTLGY TFDTASSLIN HIKL
Source	Mammalian cell
Target Names	Os08g0218700
Protein Names	Recommended name: Uncharacterized protein Os08g0218700/LOC_Os08g12160 Alternative name(s): Unknown protein AN01 from 2D-PAGE of anther
Expression Region	31-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.