



Recombinant *Oryza sativa* subsp. *japonica* Germin-like protein 8-7 (GER6)

| | |
|--------------------------|---|
| Product Code | CSB-YP754519OFG |
| Abbreviation | GER6 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6YZZZ |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | SDPSPLQ DFCVADKHSP VLVNGFACLD PKYVNADHFF KAAMLDTPRK TNKVGSNVTL INVMQIPGLN TLGISIARID YAPLGENPPH THPRATEILT VLEGTLYVGF VTSNPNTLF SKVLNKGDFV VFPEGLIHFFQ FNPNP HQPAV ALAALSSQNP GAITIANAVF GSKPPISDDI LAKAFQVEKG TIDWLQAQFW ENNHY |
| Source | Yeast |
| Target Names | GER6 |
| Protein Names | Recommended name: Germin-like protein 8-7 Alternative name(s): Germin-like protein 6 Short name= OsGER6 |
| Expression Region | 24-225 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |