



Recombinant *Oryza sativa* subsp. *japonica* COBRA-like protein 1 (BC1L6)

| | |
|--------------------------|--|
| Product Code | CSB-EP754523OFG |
| Abbreviation | BC1L6 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q6Z4G8 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | YD PLDPNGNITI KWDVLQWTPD GYVAVVSLYN YQQYRHIQSP GWKLGWVWAK KEIIWAMNGG QATEQGDCSK FKSNIHCCK KDPEIVDLLP GTPYNMQIAN CCKGGVLNSW AQDPANAIAS FQVSVGQAGT TNKTVRVP FTLKSPGPGY TCGSAKVVRP TKFFSQDGRR TTQAHMTWNV TCTYSQIVAQ RSPTCCVSL S FYNDTIVNC PTCSCGCQNN KPGSCVEGNS PYLASVVNTH NKDSLTPLVQ CTSHMCPIRV HWHVKVNYKE YWRVKITVTN FNYRMNYSQW NLVTQHPSFD NLTTIFSFNY KSLNPYGVIN DTAMLWGIKY YNDLLMTAGP DGNVQSELLF KKDKPSFTFE KGWAFPRRVY FNGDNCVMPP PDAYPWL PNA |
| Source | E.coli |
| Target Names | BC1L6 |
| Protein Names | Recommended name: COBRA-like protein 1 Alternative name(s): Protein BRITTLE CULM1-like 6 |
| Expression Region | 29-420 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |