



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 25 (Os02g0685600, LOC_Os02g46080)

Product Code	CSB-BP754533OFG
Abbreviation	Os02g0685600, LOC_Os02g46080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZHC8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MFSWLLRIAS ACLGPARRYA RTRKDEDGGD NGGGVADGLL WSRDLGRHAA GEFSFAVVQA NEALEDHSQV ETGSAATFVG VYDGHGGADA ARFISDHLFA HLIRLARESE TVSEEVVRGA FSATEEGFLT LVRRTQFLKP MIAAVGSCCL VGIIWRGVLY VANLGDSRAV VGYLGRNKI TAEQITRDHN ACKEEVRQEL ISRHPDDSQI VVLKHGVWRI KGIIQVSRTI GDAYLKRREF ALDPSITRFR LSEPLRRPVL TAEPSICTRV LSLQDQFVIF ASDGLWEHLT NQQAVIDIVYK NPRAGIAKRL VNTALKEAAR KREMRFVDLK KVEKGVRRFF HDDITVVVVY IDHELLQEKN VSVPELSVRG FVDSVGPSRI SGFDAIS
Source	Baculovirus
Target Names	Os02g0685600
Protein Names	Recommended name: Probable protein phosphatase 2C 25 Short name= OsPP2C25 EC= 3.1.3.16
Expression Region	1-387
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.